KaVo. Dental Excellence.

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Ref.: 1007 2981 edição 01

Pantone 1525 C

1 APPLICATION.....	06
2 PRINCIPLE OF OPERATION.....	06
3 PRODUCT PRESENTATION.....	06
4 PERFORMANCE OF THE DEVICE.....	06
5 PRECAUTIONS, RESTRICTIONS AND SPECIAL WARNINGS.....	08
5.1 Materials to be used with this device.....	08
5.2 LED replacement.....	08
5.3 Care during transportation.....	08
5.4 Care during use.....	08
5.5 Disposal of equipment.....	08
5.6 Biocompatibility.....	08
5.7 Chemical, physical and biological properties.....	08
5.8 Parts, accessories and materials to be used with the product.....	08
5.9 Intentional radiation.....	08
5.10 Unintentional radiation.....	08
5.11 Alterations in the operation of the LED light curing system.....	08
5.12 Battery.....	08
6 SPECIFICATIONS.....	10
7 TRANSPORTATION AND STORAGE.....	12
8 DESCRIPTION OF COMPONENT PARTS.....	14
9. SYMBOLS.....	14
10 OPERATION.....	16
10.1 Attaching the tip.....	16
10.2 Attaching the eye protector.....	16
10.3 First change.....	16
10.4 Checking the charge.....	16
10.5 Battery.....	16
10.6 Auto power off.....	16
10.7 Polymerization cycles.....	16
11. AUDIBLE WARNINGS.....	18
12. IMPORTANT NOTES.....	18
12.1 Statutory provisions.....	18
12.2 Normative classification of the product.....	18
13 CARE AND CALIBRATION.....	20
13.1 Fiber optic tip.....	20
13.2 Eye protection.....	20
13.3 Hand piece.....	20
14 ELIMINATING FAILURES.....	20
15 INSTALLATION.....	22
15.1 Place of installation.....	22
16 MAINTENANCE.....	22
16.1 Important notes.....	22
WARRANTY TERMS.....	24

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## Warranty Terms

§1. The product is guaranteed against manufacturing defects within 12 months (including the legal warranty time) as of the delivery date of the product, it would be necessary to present the Invoice and a filled Warranty Certificate.

§ 2. In case of a special extension of the warranty deadline, it does not cover the rubber components, for example, ring, diaphragms, seals and hoses.

§ 3. The LEDs of membrana are not covered by warranty.

§ 4. The guarantee is limited to repairing or replacing parts with manufacturing defects, proved by the manufacturer or Authorized Service, not including the repair of defects originating from:

- Disregard of instructions for use and maintenance;
- Falls or shock;
- Improper storage;
- Action of natural agents;
- Technical assistance carried out by persons not authorized by KaVo;
- Damage to painted parts and plastic parts (painted or not) caused by misuse, dental materials as acrylic, gloves, etc. and sharp or point end objects;
- Connection of the device to power supply other than that specified for operation thereof;
- Lack of grounding or inefficient grounding;
- Use of accessories, spare parts and materials not specified by KaVo.

§ 5. The warranty shall automatically become void:

- At its expiration date;
- For changes of the equipment or inappropriate use of accessories by the user;
- By changing the Certificate of Warranty;
- For technical assistance carried out by persons not authorized by KaVo.

§ 6. Repair or replacement of parts during warranty period does not extend the period of original validity.

§ 7. Parts replaced due to warranty, become property of KaVo of Brazil Ind. e Com Ltda.

§ 8. Buyer will bear charges due to (regardless of the product warranty):

- Installation of the product (which will be informed at sale of the product);
- The transport of parts or of the product;
- Service personnel rates at calls for assistance as of the seventh month during warranty period;
- From offset expenses of personnel involved in the care of the call to installation or technical assistance, based on the value of distance traveled (in km) to be informed in the Customer Support Network KaVo;

-The labor, materials, parts and necessary adequations for site preparation for installation of product, ie, power supply, air, water, drainage, suction, grounding and so on.

§ 09. The KaVo of Brazil Ind. e Com. does not authorize any person or entity to take on its behalf, any responsibility related to its products than those explained here.

§ 10. The buyer must sign the Certificate of Warrantee at the time of purchase.

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## ENGLISH

**1 APPLICATION**

The LED light curing system is a device used in the photo polymerization of composite resins, copolymers and glass ionomer materials. To ensure a safe polymerization, these materials must use camphorquinone as the photo initiator. Other dental materials that employ other types of photo initiators, whose absorption spectrum is outside the range of 420-480nm, are not compatible with this photo polymerizer.

This product should only be used by a trained user, and the user should observe the current regulations for the prevention of accidents and safety at work and these instructions. It is the obligation of the Users to use work instruments that are in perfect condition, observe the correct use and protect himself, the patients and others from any possible danger.

The LED light curing system is also indicated for activation/acceleration in the treatment of dental whitening.

**2 PRINCIPLE OF OPERATION**

This is based on a blue LED light- with a wavelength between 420-480 nm, used to polymerize photosensitive substances that use camphorquinone as the photo initiator.

Camphorquinone, is a common photo initiator used in dentistry materials, it has a maximum absorption within this range (465 nm). When the camphorquinone is exposed to light, in the presence of co-initiators such as amines, for example, radicals are formed which, in turn, initiate the polymerization process.

The device is composed of a source of generating a pure blue light (LED - light emitting diode), located in a hand piece (b) and controlled by an electronic circuit. The device is composed of a light source that generates a pure blue light (LED - light emitting diode), located in a hand piece (b) and controlled by an electronic circuit. The power supply uses a large capacity LITHIUM-ION battery. The hand piece (b) is portable and has no electrical wiring. The charger (c) requires electrical energy for it to operate, it should be connected to a 110 - 220 VAC, 50 - 60 Hz source.

As opposed to halogen lamps, the light emitting diode (LED) does not produce visible light by means of heating metallic filaments. As excessive heat is not produced, an internal fan is not required.

Compared with conventional lamps, the light produced by the LED has a narrower spectral distribution. This is the main difference between the light produced by the LED and that produced by halogen lamps. The LED with appropriate power only produces the wavelength range desired.

**3 PRODUCT PRESENTATION**

- a. A User Manual with the Warranty Certificate
- b. Hand piece
- c. Charger
- d. Eye protector
- e. Fiber Optic Tip

**4 PERFORMANCE OF THE DEVICE**

So that the use and handling of the LED light curing system does not cause risks to the user and the people who are in the place where it is being used, the NBR IEC 60.601-1 standard was applied in the project and the production of the device.

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
## 15 INSTALLATION

## 15.1 Installation location

⚠ It is the buyer's responsibility to provide suitable wiring (including an appropriate ground) from the connection point to the KaVo product, with the services of duly trained personnel.

If the buyer chooses to use the Authorized Service KaVo for the installation, the buyer will be responsible setting up an account for the installation and other expenses. Charges for travel and accommodation of personnel involved in the installation call, outside the urban perimeter of the city where the authorized service is located, will also be paid for by the buyer.

If the technician detects a lack of a ground or an inefficient ground, the installation of the product ceased. Once these instructions are received, the customer can contact KaVo authorized service for more details.

## 16 MAINTENANCE

## 16.1 Important notes

⚠ The maintenance should be carried out only by authorized KaVo technicians, being correctly authorized and identified.

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
d. Protetor ocular

## ENGLISH

## 5 PRECAUTIONS, RESTRICTIONS, WARNINGS AND SPECIAL CARE

**5.1 Materials to be used together with the equipment**  
Dental materials that do not contain camphoroquinone based photo initiators with maximum absorption of the wavelength range between 440 and 480 nm are not compatible with the LED light curing system.

**5.2 Replacing the LED**

In spite of having a long useful life, the LED may burn out. In order to offer the greatest safety to the equipment the LED should only be replaced by an authorized KaVo technician.

**5.3 Care in transportation**

The unit should not be dropped or bumped, because of the risk of breaking components.

**5.4 Care in use**

The device must not be dropped or bumped, because of the risk of breaking components. Special care must be given to the tip.

**5.5 Disposal of the equipment**

At the end of the useful life of the product, it must be disposed of in the appropriate place (according to local legislation).

**5.6 Biocompatibility**

The materials used in the making of the hand piece are considered harmless when in contact with the skin for less than 24 hours.

**5.7 Chemical, physical and biological properties**

The LED light curing system is non invasive and does not have a specification of chemical or biological interaction with live tissue. The system does not generate heat so there are no risks of flammability and it does not emit or lose substances.

**5.8 Parts, accessories and materials to be used with this product**

There is no prescribed use of other parts, accessories or materials besides those mentioned in these instructions. ⚠ KaVo will not assume any responsibility for damage caused by using accessories, replacement parts and materials not specified by the company.

**5.9 About intentional radiation**

The LED light curing system does not have dangerous levels of radiation, as the radiation power is distributed over an area, providing intentional radiation.

**5.10 About unintentional radiation**

The LED light curing system has a tip, which directs the light to the desired area, there is no reflected, diffused or escaped radiation.

**5.11 Alterations in the operation of the LED light curing system**

If there is any changes when using this device, stop using it. Turn it off and contact the authorized KaVo service network.

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## 13 CARE AND CALIBRATION

### 13.1 The fiber optic tip

- The fiber optical tip should be cleaned with water and soap or another mild soapy solution.
- Disinfection can be carried out using a disinfecting solution that does not have an acidic, ammonia or chlorine base.
- Sterilization can be done in a steam autoclave at 135 ° C.
- Do not use any instrument or abrasives on the sides of the tip, as they may scratch it, causing a loss in the light output.
- If the tip is damaged in any way, such as by falling, it should be replaced by a new one, as this can affect the performance of the product.

### 13.2 Eyecup

- The eye protector should be cleaned with neutral soap and water. Do not use ethyl alcohol.
- Disinfection and sterilization are not applicable.

### 13.3 Hand piece

- The hand piece should never be submerged in any liquid;
- The hand piece must be cleaned superficially with a neutral detergent solution, which has no ammonia or acid base.
- The device should have its surfaces disinfected with a neutral disinfectant solution, that is not acid or ammonia based.
- The hand piece should not be sterilized.

#### ⚠ ATTENTION

KaVo is not responsible for any damage caused by products/levels of concentration not recommended. The LED light curing system unit was designed to undergo maintenance work when needed there are no risks of it going out of calibration by the lack of maintenance.

## 14 ELIMINATING FAILURES

The corrective measures described here, can be carried out by the user. If these solutions are not sufficient to solve the problem, request the services of an authorized KaVo technician, in order to eliminate the failure. By means of an agreement with the user, KaVo will retain the circuit diagrams, lists of component, descriptions, calibration and verification instructions and other information necessary for the user's technical staff, for repairing the parts of the equipment that are considered to be repairable.

Failures	Probable cause	Solution
There is no light emission or the light level is low when turned on	Empty or partially discharged battery	Recharge the battery
The equipment is not polymerizing the resins	Resin is not suitable for the wavelength range of the LED light curing system	Acquire the appropriate resin for the wavelength of the LED light curing system i.e. which contains photo initiators with camphoroquinone

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## ENGLISH

## 11 AUDIBLE WARNINGS

- 1 beep every 10 seconds** = time indicator  
**4 consecutive beeps** = low battery indicator

## 12 IMPORTANT NOTES

## ⚠ ATTENTION

- Always follow the manufacturers' instructions on cure times for the various materials used, remember that it is better to increase the cure times than use shorter times than those recommended, because by exceeding the recommended cure time will not harm the restoration.
- Doubling the cure period does not necessarily mean that you are doubling its strength.
- The end of the optical probe should be as close as possible to the surface being cured (avoid physical contact between with the parts), note that the flat surface should be parallel to the surface being treated.
- The LED light curing system may be used for the activation/acceleration of dental whitening.

## ⚠ ATTENTION

Only qualified KaVo trained Authorized Technical Network technicians are authorized to install and provide technical assistance for KaVo products.

KaVo will not accept any responsibility for damage caused by:

- External-effects (poor quality of source supplies or the improper installation of these).
- Use of information provided by third parties.
- Repairs performed by personnel unauthorized by KaVo.
- Use of accessories, spare parts and materials not specified by KaVo Brazil.

## 12.1 Statutory provisions

It is indispensable to apply and comply with the general directives/guidelines relevant to medical products and the laws and regulations, such as the techniques for the start up and operation of the KaVo product according to the intended purpose of the product.

## 12.2 Normative classification of the product

-In accordance with the type of protection against electrical shocks	Classe II equipment
In accordance with the grade of protection against electrical shocks	Applied part of BF
In accordance with the grade of protection against dangerous water penetration	IPX0
In accordance with the grade of application safety in the presence of an inflammable anesthetic mixture	Equipment not suitable for use in the presence of an inflammable anesthetic mixture with air, oxygen or nitrous oxide
In accordance with the operational mode	Continuous operation with intermittent load. The time of operation is between 60 seconds and the rest period is 150 seconds.

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


## ENGLISH

## 7 TRANSPORT AND STORAGE

Observe the instructions printed on the outside of the packaging that are applicable to the transport and storage.

	Transport and store the packaging respecting the indicated temperature
	Transport and store the packaging with the side indicated by the arrow pointing upwards
	Protect against blows and knocks
	Protect against from moisture
	Maximum weight allowed on the packaging
	Humidity limit
	Atmospheric pressure limit

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
## 10 OPERATION

## 10.1 Fitting the tip

Fit the tip into the socket located at the end of the hand piece.

## 10.2 Fitting the eyecup

Remove the tip from the hand piece. Assemble the eye protector directly onto the body of the tip, without affecting the correct fitting of the hand piece.

## 10.3 First charge

Connect the power cord from the charger to the outlet. A light on the charger will appear - power (red) on the left and the charging indicator (green on the right). The hand piece should be placed in the charger. During charging, the power light (red) is lit and the battery control lamp (green center) will blink indicating that the device is charging.

As the battery increases its load, the charge indicator lamp (green on the right) will light up gradually until the battery control lamp (green center) stops flashing. When the charge indicator lamp (green right) is fully lit, the charge is complete. To ensure the highest level of performance of the apparatus before use, the battery must be thoroughly charged. On average, the first recharge should take from 8 to 10 hours.

When the load is complete, remove the hand piece from the charger and unplug the power cord from the network.

## 10.4 Check the charge level

Connect the power cord from the charger into the outlet. The charger will turn the power light on (red) left and the charging indicator (green on the right). The hand piece should be placed in the charger to see the intensity of the charge indicator lamp (green on the right).

If the lamp is completely off, this indicates that the product needs to be charged.

## 10.5 Battery warning

When the User selects a cure cycle, the LED membrane that indicates the cycle will remain lit. When the LED is blinking and the device sounds 4 beeps, the battery needs to be recharged. If this occurs during a procedure, increase the cure period and recharge the battery as soon as possible.

## 10.6 Auto shutdown

In order to avoid wasting battery power, the device has an auto off timer, which operates when the product has been idle for 30 seconds. To turn it on again, press the switch ①.

## 10.7 Cure cycles (polymerization)

a) 10 seconds Cycle: Select the time of 10s with the ① switch. The equipment will start the light emission with a power of 1100 mW / cm<sup>2</sup> for three seconds. Then the power will gradually reduce for three seconds to reach a power of 600 mW / cm<sup>2</sup>, remaining constant until the end of the cycle. The emission of the light can be stopped at any time by pressing the switch ①.

b) 20 seconds cycle: Select the time of 20s with the ① switch. The equipment will start the emission of light with a preheating time of seven seconds, gradually increasing the power to reach a peak of 1100 mW / cm<sup>2</sup>. It will remain like this for three seconds, and then gradually reduce for another three seconds, until the reaching a power of 600 mW / cm<sup>2</sup>, remaining constant until the end of the cycle. The light emission can be stopped at any time by pressing the switch ①.

c) 40 seconds cycle: Select the time of 40s with the ① switch. The equipment will start the emission of light with a preheating time of seven seconds, gradually increasing the power to reach a peak of 1100 mW / cm<sup>2</sup>. It will remain like this for three seconds, and then gradually reduce for another three seconds, until the reaching a power of 600 mW / cm<sup>2</sup>, remaining constant until the end of the cycle. The light emission can be stopped at any time by pressing the switch ①.

d) 60 seconds cycle: Select the time of 60s with the ① switch. The equipment will start the emission of light with a preheating time of seven seconds, gradually increasing the power to reach a peak of 1100 mW / cm<sup>2</sup>. It will remain like this for three seconds, and then gradually reduce for another three seconds, until the reaching a power of 600 mW / cm<sup>2</sup>, remaining constant until the end of the cycle. The light emission can be stopped at any time by pressing the switch ①.

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